

## Pharmacomechanical thrombolysis treatment of DVT and preservation of valve competency

Venous thromboembolism (VTE), which includes both deep vein thrombosis (DVT) and/or pulmonary embolism (PE), affects approximately 5% of the population at some point in their lives<sup>1</sup> and is estimated to occur in 300,000 – 600,000 persons in the United States each year.<sup>2</sup> Most patients undergo systemic anticoagulation and compression stocking therapy, a treatment strategy intended to minimize propagation of the clot and prevent pulmonary embolization.<sup>1</sup> The thrombus is not removed, but rather relies upon the endogenous fibrinolytic system to eventually resolve the thrombus and recanalize the vein.<sup>3</sup> However when treated with anticoagulation alone, it is estimated that in only 12-47% of cases is there complete lysis the offending clot within 6 months,<sup>4,5</sup> increasing the risk of recurrent DVT and/or post thrombotic syndrome (PTS),<sup>5,8</sup> a debilitating late sequelae of DVT.

Post thrombotic syndrome (PTS) occurs in nearly 30-50% of patients following an episode of DVT<sup>8,9</sup> and is characterized by chronic leg pain and swelling that limits activity, the ability to work, and potentially results in leg ulcers.<sup>9</sup> PTS is thought to be caused by venous hypertension and impaired venous return due to the persistent (residual) venous obstruction.<sup>9,10</sup> Valvular reflux also frequently occurs after DVT, likely via thrombus-induced activation of inflammation, fibrous scarring by acute and resolving thrombus, or by venous dilation distal to the obstructed segment.<sup>8</sup> Risk factors of PTS include symptom severity at 1-month post DVT diagnosis, thrombosis of the common femoral or iliac vein, higher body mass index (BMI), previous ipsilateral venous thrombosis, older age, and female gender.<sup>11</sup> The occurrence of PTS has been shown to be directly related to the degree of residual clot; conversely, in cases where when there is a high degree of clot lysis, rethrombosis is avoided and patients have a minimal risk of PTS.<sup>3,12</sup> Furthermore, it has been recently shown that the use of compression stockings do not prevent the occurrence of PTS after DVT nor do they influence the severity of PTS or the rate of recurrent VTE.<sup>13</sup>

Including early clot removal in the treatment strategy for management of DVT has been shown to be associated with a significant reduction in the occurrence of persistent venous obstruction, recurrent thrombosis, PTS, valvular incompetence (reflux) and improved quality of life (QOL),<sup>1,4,14</sup> presumably by rapidly restoring venous patency and preserving or limiting damage to venous values.<sup>15</sup>

Due to the significant impact of both PTS and recurrent thrombosis on QOL and productivity as well as the cost associated with managing these conditions,<sup>11,14</sup> the Society for Vascular Surgery (SVS) has recently updated clinical practice guidelines suggesting early removal of thrombus in patients most likely to benefit from the procedure.<sup>16</sup> As defined by the SVS, early thrombus removal is indicated for patients with:

- A first episode of acute iliofemoral DVT
- Symptoms of <14 days in duration
- Low bleeding risk
- Ambulatory with good functional capacity and an acceptable life expectancy, who would benefit from a reduced long-term risk of PTS and DVT recurrence

First-line therapies for early thrombus removal include catheter-directed pharmacologic thrombolysis (CDT) and pharmacomechanical thrombolysis (PMT). Catheter-directed thrombolysis involves infusion of a thrombolytic agent in the region of the venous thrombus through a percutaneously placed catheter. Catheter-directed thrombolysis is limited by the risk of bleeding (5-11%),<sup>6</sup> prolonged infusion time required for the thrombolysis procedure (average 48-53 hours),<sup>17</sup> and the cost associated with the frequent requirement of patient monitoring in the ICU. Pharmacomechanical thrombolysis (PMT) was developed to address these limitations through the adjunctive use of mechanical devices, involving rotational, rheolytic or ultrasound technologies, to deliver the thrombolytic agent as well as produce some combination of thrombus fragmentation, distribution of thrombolytic drug through the thrombus and/or thrombus aspiration.<sup>18</sup>

Comparisons of CDT with PMT have shown that PMT is associated with<sup>20,21</sup>:

- Comparable rates of treatment success and risk of complications
- Shorter treatment times
- Lower doses of thrombolytic drug required
- Reduced ICU and total hospital stay and hospital costs

Recognizing the benefits of PMT over CDT, the SVS now suggests the use of PMT over catheter-directed thrombolysis alone if expertise and resources are available.<sup>16</sup>

Despite these advantages, concerns remain that mechanical device may damage vein valves as compared to CDT. The effect of PMT on venous valve function as compared with conventional CDT therapy has been recently evaluated by Vogel and colleagues.<sup>14</sup> Post procedural phlebograms were reviewed to determine residual obstruction; valve function (reflux) was evaluated using duplex ultrasonography at follow-up. A total of 20 limbs were evaluated that were treated with CDT alone, 49 were treated with PMT (n=35 with the Trellis catheter, n =14 with AngioJet™ System). Similar rates of success were achieved with both CDT and PMT. At a mean length of follow-

**Table 1: Summary of literature evaluating venous valve function following treatment for DVT**

Author	N	DVT Therapy	Method for Determination of Venous Valve Closure Time	Length of Follow-up (months)	% with Normal Valve Function	p
Vogel, 2012	20 49	CDT PMT	Duplex Ultrasonography	44	CDT: 35% PMT: 47%	.42
Arko, 2007	30	PMT	Duplex Ultrasonography	12	88%	--
Elsharawy, 2002	17 18	Anticoagulation CDT	Photoplethysmography	6	Anticoag: 59% CDT: 89%	0.04
VanHaarst, 1996	24	Anticoagulation	Duplex Ultrasonography	33 85	52% 40%	--

up of 44 months, there was no difference in the rate of reflux demonstrated between patients treated with CDT (65%) and PMT (53%) ( $p=.42$ ). Limbs treated with the Trellis catheter exhibited a reflux rate of 57% vs. 43% in limbs treated with AngioJet System ( $p=.53$ ). Although the sample size was insufficient to demonstrate statistical significance, the risk of valve reflux following treatment with PMT was shown to be at least equivalent, and possibly slightly lower, than that which occurs following CDT.

Studies have reported a wide variation in rates of reflux after the use of CDT or anticoagulation alone, depending on the method of assessment and length of follow-up, as shown in Table 1. The results of the Vogel series are consistent with other reports evaluating valvular function following treatment of DVT.

The AngioJet System is a rheolytic PMT technology that is used extensively to treat DVT, acute limb ischemia and thrombosis of hemodialysis access grafts. Data collected on 329 lower extremity DVT patients as a part of the PEARL Registry, a study following AngioJet patients funded by Boston Scientific has shown that:

- Thrombectomy with Power Pulse™ or rapid lysis was utilized in 87% of the patients.
- An average reduction of 96% in thrombus immediately post AngioJet PMT therapy.
- Rates for freedom from rethrombosis at 6 and 12 months were 88% and 84%, respectively.
- Clinically and statistically significant improvements in both physical and mental components of QOL measures from baseline to 3, 6, and 12 months.
- 75% of AngioJet cases were completed in  $\leq 24$  hours.
- Procedure-related complications included bleeding requiring transfusion (N=1), hematoma at access site (2) arrhythmia during the AngioJet procedure (3), Acute renal failure (1) or Transient renal failure (1), and pulmonary embolism (1).

As shown in Table 2 below, the lower extremity outcomes for AngioJet PMT therapy as reported in the PEARL registry aligned favorably with the CDT experience as reported by the Venous Registry<sup>17</sup> with considerably less lytic drip times when CDT is used in combination with AngioJet therapy. In addition, the PEARL registry reported greater thrombus removal both acute and chronically with fewer bleeding complications with AngioJet Thrombectomy use in the treatment protocol.

**Table 2: Comparison of PEARL Registry and Venous Registry patient populations and outcomes in the treatment of lower extremity DVT**

	Venous Registry <sup>17</sup>	PEARL Registry
Number of Patients	287	329
# of sites	63	32
Prior DVT	31%	40%
Concomitant stent placement	33%	35%
Primary access	Popliteal	Popliteal
Mean age	47.5 years	52 years
% Male	48%	57%
Primary treatment	CDT	PMT with or without CDT
Primary lytic	Urokinase	TPA
CDT drip times (mean)	48 hours	17 hours (when CDT was utilized)
Procedure time (median)	NA	CDT: 41 hours CDT+PMT: 22 hours PMT: 2 hours
Acute: % thrombus removal	86%	97%
Chronic: % thrombus removal	68%	95%
Primary Patency	6 month: 65% 12 month: 60%	NA
Freedom from Rethrombosis	NA	6 month: 88% 12 month: 84%
Bleeding complications	11% major 16% minor	5% (major + minor)

## Conclusions

While definitive, randomized trials evaluating the late outcomes for treatment of DVT using anticoagulation alone vs. aggressive early intervention with PMT are ongoing,<sup>24</sup> the consistent results reported by an accumulating, large body of research comprised of smaller studies and registries supports the benefit of early thrombus removal using PMT:

- Catheter based thrombus removal for DVT treatment has been shown to improve patient outcomes.
- Advantages of PMT treatment of DVT are not compromised by adverse effects on vein valve function.
- Studies support PMT as a safe method for the treatment of symptomatic acute DVT with excellent rates of patency reported.
- Rates of recurrent venous thromboembolism have been shown to be lower in treatment groups receiving percutaneous venous intervention.

## References

- 1 Casey ET, Murad MH, Zumaeta-Garcia M, et al. Treatment of acute iliofemoral deep vein thrombosis. *J Vasc Surg.* 2012; 55:1463-73.
- 2 Beckman MG, Hooper WC, Critchley SE, Ortel TL. Venous thromboembolism: a public health concern. *Am J Prev Med.* 2010; 38:S495-501
- 3 Comerota AJ, Grewal N, Martinez JT, et al. Postthrombotic morbidity correlates with residual thrombus following catheter-directed thrombolysis for iliofemoral deep vein thrombosis. *J Vasc Surg.* 2012; 55:768-73.
- 4 Elsharawy M, Elzayat E. Early results of thrombolysis vs anticoagulation in iliofemoral venous thrombosis. A randomised clinical trial. *Eur J Vasc Endovasc Surg.* 2002; 24:209-14.
- 5 Enden T, Haig Y, Kløw NE, et al.; CaVenT Study Group. Long-term outcome after additional catheter-directed thrombolysis versus standard treatment for acute iliofemoral deep vein thrombosis (the CaVenT study): a randomized controlled trial. *Lancet.* 2012; 379:31-8.
- 6 Comerota AJ. Thrombolysis for deep venous thrombosis. *J Vasc Surg.* 2012; 55:607-11.
- 7 Hull RD, Marder VJ, Mah AF, Biel RK, Brant RF. Quantitative assessment of thrombus burden predicts the outcome of treatment for venous thrombosis: a systematic review. *Am J Med.* 2005;118(5):456-64.
- 8 Kahn SR. The post-thrombotic syndrome. *Hematology Am Soc Hematol Educ Program.* 2010; 2010:216-20.
- 9 Prandoni P, Kahn SR. Post-thrombotic syndrome: prevalence, prognostication and need for progress. *Br J Haematol.* 2009;145:286-95.
- 10 Vedantham S. Valvular dysfunction and venous obstruction in the post-thrombotic syndrome. *Thromb Res.* 2009;123 Suppl 4:S62-5.
- 11 Kahn SR, Shrier I, Julian JA, et al. Determinants and time course of the postthrombotic syndrome after acute deep venous thrombosis. *Ann Intern Med.* 2008; 149:698-707.
- 12 Grewal NK, Martinez JT, Andrews L, Comerota AJ. Quantity of clot lysed after catheter-directed thrombolysis for iliofemoral deep venous thrombosis correlates with postthrombotic morbidity. *J Vasc Surg.* 2010; 51:1209-14.
- 13 Kahn, SR, Shapiro S, Wells PH, et al. A multicenter randomized placebo controlled trial of compression stockings to prevent the post-thrombotic syndrome after proximal deep venous thrombosis. *Proceedings of the American Society of Hematology;* December, 2012; Atlanta, GA. Abstract 393.
- 14 Ashrani AA, Silverstein MD, Rooke TW et al. Impact of venous thromboembolism, venous stasis syndrome, venous outflow obstruction and venous valvular incompetence on quality of life and activities of daily living: a nested case-control study. *Vasc Med.* 2010; 15:387-97.
- 15 Comerota AJ, Paolini D. Treatment of acute iliofemoral deep venous thrombosis: a strategy of thrombus removal. *Eur J Vasc Endovasc Surg.* 2007; 33:351-60.
- 16 Meissner MH, Gloviczki P, Comerota AJ, et al; Society for Vascular Surgery; American Venous Forum. Early thrombus removal strategies for acute deep venous thrombosis: clinical practice guidelines of the Society for Vascular Surgery and the American Venous Forum. *J Vasc Surg.* 2012; 55:1449-62.
- 17 Mewissen MW, Seabrook GR, Meissner MH, et al. Catheter-directed thrombolysis for lower extremity deep venous thrombosis: report of a national multicenter registry. *Radiology.* 1999; 211:39-49.
- 18 Karthikesalingam A, Young EL, Hinchliffe RJ, et al. A systematic review of percutaneous mechanical thrombectomy in the treatment of deep venous thrombosis. *Eur J Vasc Endovasc Surg.* 2011;41:554-65.
- 19 Lin PH, Zhou W, Dardik A, et al. Catheter-direct thrombolysis versus pharmacomechanical thrombectomy for treatment of symptomatic lower extremity deep venous thrombosis. *Am J Surg.* 2006; 192:782-8.
- 20 Rao AS, Konig G, Leers SA, et al. Pharmacomechanical thrombectomy for iliofemoral deep vein thrombosis: an alternative in patients with contraindications to thrombolysis. *J Vasc Surg.* 2009; 50:1092-8.
- 21 Martinez Trabal JL, Comerota AJ, et al. The quantitative benefit of isolated, segmental, pharmacomechanical thrombolysis (ISPMT) for iliofemoral venous thrombosis. *J Vasc Surg.* 2008; 48:1532-7
- 22 Vogel D, Walsh ME, Chen JT, Comerota AJ. Comparison of vein valve function following pharmacomechanical thrombolysis versus simple catheter-directed thrombolysis for iliofemoral deep vein thrombosis. *J Vasc Surg.* 2012; 56:1351-4.
- 23 van Haarst EP, Liasis N, van Ramshorst B, Moll FL. The development of valvular incompetence after deep vein thrombosis: a 7 year follow-up study with duplex scanning. *Eur J Vasc Endovasc Surg.* 1996;12:295-9.
- 24 Vedantham S, Goldhaber SZ, Kahn SR, et al. Rationale and design of the ATTRACT Study: A multicenter randomized trial to evaluate pharmacomechanical catheter-directed thrombolysis for the prevention of postthrombotic syndrome in patients with proximal deep vein thrombosis. *Am Heart J.* 2013;165:523-530.

#### ANGIOJET™ ULTRA CONSOLE

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

**INTENDED USE/INDICATIONS FOR USE:** The Console is intended for use only in conjunction with an AngioJet Ultra Thrombectomy Set. Refer to the individual Thrombectomy Set Information for Use manual for specific clinical applications. **CONTRAINDICATIONS:** Refer to the individual Thrombectomy Set Information for Use manual for specific contraindications. **WARNINGS and PRECAUTIONS:**

- Use the AngioJet Ultra Console only with an AngioJet Ultra Thrombectomy Set. This Console will not operate with a previous model pump set and catheter.
- Do not attempt to bypass any of the Console safety features.
- If the catheter is removed from the patient and/or is inoperative, the waste tubing lumen, guide catheter, and sheath should be flushed with sterile, heparinized solution to avoid thrombus formation and maintain lumen patency. Reprime the catheter by submerging the tip in sterile, heparinized solution and operating it for at least 20 seconds before reintroduction to the patient.
- Refer to the individual AngioJet Ultra Thrombectomy Set Information for Use manual for specific warnings and precautions.
- Do not move the collection bag during catheter operation as this may cause a collection bag error.
- Monitor thrombotic debris/fluid flow exiting the catheter through the waste tubing during use. If blood is not visible during console activation, the catheter may be occlusive within the vessel or the outflow lumen may be blocked.
- Ensure adequate patient anticoagulation to prevent thrombus formation in outflow lumen.
- Refer to individual Thrombectomy Set Instructions for Use manual for specific instructions regarding heparinization of the Thrombectomy Set.
- The Console contains no user-serviceable parts. Refer service to qualified personnel.
- Removal of outer covers may result in electrical shock.
- This device may cause electromagnetic interference with other devices when in use. Do not place Console near sensitive equipment when operating.
- Equipment not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Where the "Trapping Zone Hazard for Fingers" symbol is displayed on the console, there exists a risk of trapping or pinching fingers during operation and care must be exercised to avoid injury.
- Do not reposition or push the console from any point other than the handle designed for that purpose. A condition of overbalance or tipping may ensue.
- The AngioJet Ultra Console should not be used adjacent to or stacked with other equipment, and if adjacent or stacked use is necessary, the AngioJet Ultra Console should be observed to verify normal operation in the configuration in which it will be used.
- Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.
- The use of accessories and cables other than those specified, with the exception of accessories and cables sold by Bayer HealthCare as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the Ultra Console.
- MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding Electro-Magnetic Compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in the tables provided in the IFU.

**ADVERSE EVENTS:** Refer to the individual Thrombectomy Set Information for Use manual for specific observed and/or potential adverse events.

ANGIOJET is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates.

**Boston  
Scientific**

Advancing science for life™

#### Peripheral Interventions

300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2015 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-279906-AA JUL2015